

## Chapter 8 Optical Interferometry

- 8.1 The Michelson Interferometer
- 8.3 Variations of the Michelson Interferometer
- 8.4 The Fabry-Perot Interferometer
- 8.5 Fabry-Perot Transmission: The Airy Function
- 8.7 Variable-Input Frequency FP Interferometers

### 8.1 The Michelson Interferometers


➤ Amplitude-division interferometer:

Michelson Interferometer: *Make use of two beam interference*

Fabry-Perot (FP) Interferometer: *Make use of multiple beam interference*

➤ Applications : measure wavelength, measure thickness, flatness, etc.

Michelson Interferometers (using 2 Mirrors and 1 Beam splitter)



Compensate plate C: to compensate path difference between two paths

## 8.1 The Michelson Interferometers





Figure 9.24 The Michelson Interferometer. (a) Circular fringes are centered on the lens. (b) Top view of the interferometer showing the path of the light. (c) A wedge fringe pattern was distorted when the tip of a hot soldering iron was placed in one arm. Observe the interesting perceptual phenomenon whereby the region corresponding to the iron's tip appears faintly yellow. (Photo by E. H.)

➤ Dark fringes:  $2d \cos \theta = m\lambda$



$\theta$ : angular radius

➤ The center ring is dark with highest order  $m$ :  $m_{\max} = \frac{2d}{\lambda}$

➤ Angular separation :  $|\Delta\theta| = \frac{\lambda_0 \Delta m}{2d \sin \theta}$



## 8.1 Th



Q: do you think the circular fringes have uniform angular separation?

Q: When  $M_2$  is moved toward  $M_1$ , what happen to the fringe ?


1. The rings shrink toward the center with the highest-order disappearing whenever  $d$  decrease by  $\lambda/2$
2. Each remaining ring broadens as more and more fringes vanish at the center until only a few fill the whole screen, which can be explained by the angular separation of the ring.

## 8.1 The Michelson Interferometers

### Example 8-1

Fringes are observed due to monochromatic light in a Michelson interferometer. When the movable mirror is translated by 0.073 mm, a shift of 300 fringes is observed. What is the wavelength of the light? What displacement of the fringe system takes place when a flake of glass of index 1.51 and 0.005 mm thickness is placed in one arm of the interferometer? (Assume that the light beam is normal to the glass surface)

## 8.3 Variation of The Michelson Interferometers \*\*



## 8.4 The Fabry-Perot Interferometers

- Fabry-Perot interferometer : Formed by 2 parallel mirrors ( $r_1, r_2$ ), separated by a distance  $d$





Figure 9.43 Fabry-Perot etalon.

## 8.5 FP Transmission: The Airy function

See p186

$$\cos \delta = 1 - 2 \sin^2 \frac{\delta}{2}$$

- For lossless mirrors, the transmittance of FP cavity:

$$T = \frac{(1-r^2)^2}{(1+r^4) - 2r^2 \cos \delta} = \frac{1}{1 + \frac{4r^2}{(1-r^2)^2} \sin^2(\delta/2)}$$

Airy function

- Round-trip phase shift:  $\delta = 2kd$

- Parameters:


- Coefficient of finesses  $F$ :
- Half width  $\delta_{1/2}$
- Finesse  $\mathcal{F}$

- Coefficient of finesses:

$$F = \frac{4r^2}{(1-r^2)^2} \Rightarrow T = \frac{1}{1 + F \sin^2(\delta/2)}$$

$F$  is a measure of fringe contrast :

$$\text{Fringe Contrast} = \frac{T_{\max} - T_{\min}}{T_{\min}} = F$$




## 8.5 FP Transmission: The Airy function

- **Half width  $\delta_{1/2}$ :** *Half-width at half max. (HWHM)*  
**Fall width  $2\delta_{1/2}$ :** *Fall-width at half max. (FWHM)*

$$T_{\max} = 1; \quad \text{Half-max: } T = 1/2$$

$$T = \frac{1}{1 + F \sin^2 \frac{\delta}{2}} \quad \Leftrightarrow \quad F \sin^2 \frac{\delta}{2} = 1$$



**Two solutions**  $\delta_{x,y} = 2m\pi \pm 2 \sin^{-1} \frac{1}{\sqrt{F}} \quad \Leftrightarrow \quad \delta_x - \delta_y = 4 \sin^{-1} \frac{1}{\sqrt{F}} = 2\delta_{1/2} \quad \Leftrightarrow \quad \delta_{1/2} = 2 \sin^{-1} \frac{1}{\sqrt{F}}$

- **Finesse  $\mathcal{F}$ :** *Ratio of phase separation between adjacent transmittance peaks to the full width  $2\delta_{1/2}$*


$$\mathcal{F} = 2\pi / (2\delta_{1/2}) = \pi / \delta_{1/2} \quad \Leftrightarrow \quad \mathcal{F} = \frac{\pi}{2 \sin^{-1} \frac{1}{\sqrt{F}}}$$

Considering  $F$  is very large number:  $\mathcal{F} = \pi\sqrt{F} / 2 = \frac{\pi r}{1-r^2} \quad (\text{Eq. 8-29})$

## 8.5 FP Transmission: The Airy function



### Problem 8-2

Estimate the coefficient of finesse  $F$ , the finesse  $\mathcal{F}$ , and the mirror reflectivity  $r$  for a Fabry-Perot cavity with the transmittance curve shown in Fig. 8-10.



## 8.7 Variable-Input Frequency FP Interferometers

- When  $d$  is fixed, often call it as etalon, or FP etalon



- Resonance conditions (max condition):  $\delta = 2kd = 4\pi(\nu/c)nd = 2m\pi$

$$\nu_m = m \frac{c}{2nd}$$

- Free Spectral Range (FSR) of the FP interferometer:  $f_{FSR} = \Delta\nu = \frac{c}{2nd}$
- FSR: Frequency separation between adjacent transmittance peaks*

### Problem 8-13

Apply the reasoning used to calculate the finesse of a Fabry-Perot interferometer to the Michelson interferometer. Using the irradiance of Michelson fringes as a function of phase, calculate (a) the fringe separation; (b) the fringe width at half-maximum; (c) their ratio, the finesse.